



Durable Medical Orthotic Devices for mild to extreme contractures

HEELIFT

Made from latex free medical grade foam. Friction free tricot base for improved client mobility. Ensures Heelift stays in place at all time. Extra pad to control hip rotation, foot drop or provide added elevation.

Description	Code	Price
SMOOTH	DMHS	\$ 59.00
CONVOLUTED	DMHC	\$ 59.00



One size fits all. Machine washable and dryer safe.

HEELIFT TRACTION

Made from latex free medical grade foam. Comes with Velcro side traction straps. Includes removable traction bar and rope. Friction free tricot base for improved client mobility ensures Heelift stays in place at all time.

Description	Code	Price
SMOOTH	DMTB	\$ 117.00



One size fits all. Machine washable and dryer safe.

HEELIFT AFO

Smooth or Convolutated Foam with a semi-rigid polypropylene brace. Offloads heel pressure. Machine washable and autoclavable.

Description	Code	Price
SMOOTH	DMAFOS	\$ 130.00
CONVOLUTED	DMAFOC	\$ 130.00



ELBOWLIFT

Made from medical grade, latex free foam. Held in place with padded Velcro strap. Smooth non friction backing for improved client movement.

Description	Code	Price
SMOOTH	DMEL	\$ 25.00



Minimises shear. Centre holes creates no contact area for elbow. Also functions as a placement guide. Machine washable and dryer safe.



Heelift Suspension Boots are designed to provide an effective solution for prevention and assistive healing of pressure ulcers at the heel, while delivering positioning and alignment of the lower limb.

By using an inbuilt wedge of soft, firm medical grade foam the heel is ensured a no contact area while giving good support over the calf.

The Heelift fits comfortably onto the foot and is secured with Velcro® straps. The base of the Heelift is made up of a tricot low shear/friction surface which allows the boot to slide easily over the mattress surface. This reduces friction and shear, and ensures the Heelift stays in place at all times.

The Heelift comes with an extra foam piece which can be fitted to prevent footdrop or prevent hip rotation.


Foam can be cut to accommodate Achilles tendon erythema.

Machine washable and autoclavable.





SIZE	WAIST/HIP CIRC	SHORT LSO	AVE. LSO	TALL LSO
XSLIM	63.5cm - 76cm	30.5cm	33cm	35.5cm
SLIM	76cm - 89cm	30.5cm	33cm	35.5cm
REGULAR	89cm - 101.5cm	30.5cm	33cm	35.5cm
FULL	101.5cm - 117cm	30.5cm	33cm	35.5cm
XFULL	117cm - 134.5cm	30.5cm	33cm	35.5cm



SAFE SPINE TLSO

Uses: Comfortably treat functional curves of spine/back like Scoliosis, Lordosis, and Kyphosis. Supporting the spine and relieving pain from common conditions like osteoporosis and degenerative disc disease.

Features: Supporting patients with decreased or lack of muscle tone and certain brain and spinal cord injuries. Shoulder straps encourage shoulders and upper body toward corrected posture. Helps to relieve pain. Wear while walking or sitting.

CERVICAL EXTENSION

Uses: Treats forward head/upper body posture.

Features: Attaches to back of VertebrEase LSO, TLSO or Safe Spine TLSO. Comfortable, non rigid. Heat mouldable for best fit, and remould to gradually improve alignment. Use while walking or sitting.

HYPER EXTENSION

Uses: Treats hyper extended head leaning.

Features: Attaches to back of VertebrEase LSO, TLSO or Safe Spine TLSO. Heat mouldable for exact fit. Remould/reset to gradually correct alignment. Semi-rigid design is comfortable & applies prolonged low load passive stretch to reach relaxed state and re-lengthen shortened tissue, correcting deformities.

SAFE CERV LEFT/RIGHT

Uses: Treats left head leaning with or without forward leaning.

Features: Attaches to back of VertebrEase LSO, TLSO or Safe Spine TLSO. Heat mouldable for exact fit. Remould/reset to gradually correct alignment. Semi-rigid design is comfortable & applies prolonged low load passive stretch to reach relaxed state and re-lengthen shortened tissue, correcting deformities.

Size: Small = Head Circumference: Smaller than 55cm; Large = Head Circumference: Larger than 55cm

Length of Foot	Code	Price
Up to 7.5cm	MR20746-3	\$ 67.50
7.5cm - 9cm	MR20746-3.5	\$ 67.50
9cm - 10cm	MR20746-4	\$ 67.50



Infant / Toddler Safe Boot

Length of Foot	Code	Price
10cm - 13cm	MR20746-5	\$ 67.50
13cm - 15cm	MR20746-6	\$ 67.50
15cm - 17.5cm	MR20746-7	\$ 67.50
17.5cm - 20cm	MR20746-8	\$ 67.50



Paediatric Safe Boot

Code	Price
MR2074411	\$ 57.50



Safe Boot II

Size	Code	Price
SMALL	MR20734S	\$ 57.50
MEDIUM	MR20734M	\$ 57.50
LARGE	MR20734L	\$ 57.50



Neuroflex Restorative Boot

Code	Price
MRDAP500	\$ 180.00



Heel Z Float

Size	Womens Shoe Size	Mens Shoe Size	Code	Price
SMALL	5-7	5	LZ66426S	\$ 195.00
MEDIUM	8-11	6-9	LZ66426M	\$ 195.00
LARGE	12+	10+	LZ66426L	\$ 195.00




Turnbuckle Ankle



ADULT	
SIZE	FOOT LENGTH
Small	Up to 20cm
Medium	Up to 24cm
Large	Up to 28cm
X-Large	Up to 33cm

PEDIATRIC	
SIZE	FOOT LENGTH
Small	Up to 7.5cm
Medium	7.5cm - 9cm
Large	9cm - 10cm



1: Choose Hip Circumference 2: Choose Back Panel Length

63cm - 73cm	Extra Slim	35cm	Short
76cm - 89cm	Slim	40cm	Average
89cm - 101cm	Regular	45cm	Tall
101cm - 116cm	Full		

3: Take the two and add to the end of the code.

For example: Hip Circumference: 80cm = Slim
Back Panel: 45cm = Tall

Order Code: MR20728slim/tall

INFANT/TODDLER SAFE BOOT

Uses: Reduces tone/Spasticity.

Features: Provides the prolonged low load passive stretch to reach relaxed state and prevent shortening of tissue/deformity. Sewn on non skid sole prevents slipping during weight bearing and short ambulation. Design helps to maintain foot in alignment while flex addresses the tone.

PAEDIATRIC SAFE BOOT

Uses: For patients exhibiting foot drop (ankle) contracture, inversion/eversion, or internal/external hip rotation, with or without high degrees of tone.

Features: Provides the prolonged low load passive stretch to reach relaxed state and prevent shortening of tissue/deformity. Sewn on non skid sole prevents slipping during weight bearing and short ambulation. Design helps to maintain foot in alignment while flex addresses the tone.

SAFE BOOT II

Uses: Comfortably treats plantar flexion, ankle alignment, and foot drop.

Features: Support foot/ankle while sitting standing or walking short distances. Adjustable and removable foot extension encourages Toe Off action. Rotator bar in back helps maintain limb alignment. Breathable, moisture wicking padding. Non-skid sole attaches with Velcro. Flex Straps provide extra stretch for plantar fasciitis and to treat mild lost range and tone.

NEUROFLEX RESTORATIVE BOOT

Uses: Safely and comfortably treats inversion or eversion, foot/ankle contractures, foot drop and plantar flexion.

Features: Use whilst lying or sitting. Not suitable for walking. Suspends the heel. Heat mould to treat almost any foot/ankle deformity. Remould to gradually bring joint toward and to correction. Semi-rigid design with flex properties allows the joints to reach relaxation then continues the slow gentle stretch.

HEEL Z FLOAT

Uses: Air cushion suspends heel for reduced pressure.

Features: Allows for total suspension of the heel while maximizing immersion and pressure redistribution in other vulnerable areas. A single static air chamber utilizes air displacement to normalize pressure over the entire contact area for overall lower limb pressure relief and total heel suspension. Thermo-material maintains a comfortable temperature range to help prevent sweating and skin breakdown.

TURNBUCKLE ANKLE

Enables the clinician to apply progressive low-load torque at the ankle joint to control muscle imbalance and reduce contracture.

Code	Price
MR20728	\$ 140.00 - \$ 168.00



Safe Spine TLSO

Size	Code	Price
PAEDIATRIC	MR20748C	\$ 143.00
ADULT	MR20728C	\$ 143.00



Cervical Extension

Code	Price
MR20728H	\$ 143.00



Hyper Extension

Size	Left	Right	Price
SMALL	MR20728SCS-L	MR20728SCS-R	\$ 190.00
LARGE	MR20728SCL-L	MR20728SCL-R	\$ 190.00



Safe Cerv
Left or Right



SIZE	WIDTH OF MCP	WRIST TO FINGER
Infant	4cm-5cm	Up to 6cm
Toddler Jr.	4cm-5cm	Up to 7.5cm
Toddler	5cm—5.5cm	Up to 8.5cm
Child	5.5cm-6.25cm	Up to 10cm
Youth	6cm—7.5cm	Up to 12.5cm
X-Small	7cm - 7.5cm	Up to 13cm
Small	Up to 7.5cm	Up to 16.5cm
Medium	Up to 9cm	Up to 18cm
Large	Up to 9.5cm	Up to 20cm

RESTORATIVE HAND

Uses: Effectively and comfortably treat any wrist/hand/finger lost range and tension from fresh stroke, brain or spinal cord injuries to moderate deformities and pain.
Features: Flex design encourages joints to relax. Can heat mould, then remould to gradually correct joints to normal alignment. Provides prolonged low load passive stretch to re-lengthen shortened tissue. Kydex base extends to end of thumb. Prevents deformities if treated in early stages.

HYPER HAND PAEDIATRIC

Uses: Treats from the tightest closed palm, to moderate flexion.
Features: Soft silky padding and straps. Plastic slit between each digit can be moulded per finger.

NEUROFLEX FLEX HAND

Uses: Treat from the tightest closed palm, to moderate flexion.
Features: Efficiently works to relax joint tension, even from the wrist. Graduated sizes of cones gradually open hand/fingers/thumb.

THUMB EASE

Uses: Ideal for severe thumb deformities.
Features: Gradually add more filler to pocket under thumb to slowly bring it to better alignment. Mould Kydex base to apply only low load passive stretch, then remould to gradually take fingers/wrist from deformity toward correction. Comfortable. Completely covered with washable, breathable, moisture wicking padding with removable finger separators.

THUMB SPICA

Uses: Treats injuries like strains, sprains, & repetitive motion injuries. Excellent for patient compliance.
Features: Immobilizes thumb & wrist comfortably while leaving fingers free. Supports while joints heal. Remouldable plastic is totally covered with moisture wicking, washable padding.

COCK-UP WRIST SPLINT


Uses: Treat Carpal Tunnel Syndrome & other conditions that require stabilizing the wrist.
Features: Holds wrist in desired degree of cock-up position/shipped at 30 degrees. Easily remoulded with heat. Comfortable with soft strap through web space. Washable lining.

BEND EASE

Features: Foam over wire technology is used to comfortably and quickly address lost range of motion by easily shaping the splint to conform to the patient. Padded with Outlast thermal regulating material for added comfort and protection.

TURNBUCKLE FUNCTIONAL POSITION SPLINT

Permits application of incremental low-load force to the wrist joint to gradually reduce contracture while hand is aligned in a functional position



Size	Code	Price
PAEDIATRIC	MR20745C	\$ 180.00
X SMALL	MR20725XS	\$ 180.00
SMALL	MR20725S	\$ 180.00
MEDIUM	MR20725M	\$ 139.50
LARGE	MR20725L	\$ 139.50



Neuroflex Restorative Knee

Code	Price
MR20725SP	\$ 34.00



Neuroflex Knee Separator

Size	Length	Circ	Code	Price
SMALL	280mm	180mm	MR20741S	\$ 41.00
MEDIUM	280mm	305mm	MR20741M	\$ 41.00
LARGE	380mm	305mm	MR20741L	\$ 41.00



Safe Straight Knee

Size	Left	Right	Price
SMALL	MR20755SL	MR20755SR	\$ 191.00
MEDIUM	MR20755ML	MR20755MR	\$ 191.00
LARGE	MR20755LL	MR20755LR	\$ 191.00



Restair Knee

Size	Kneecap Circ.	Right	Price
SMALL	33cm - 38cm	LZ66425B	\$ 278.00
MEDIUM	38cm - 43cm	LZ66425C	\$ 278.00
LARGE	43cm - 45.7cm	LZ66425D	\$ 278.00
EXTRA LARGE	45.7cm - 51cm	LZ66425E	\$ 278.00



Turnbuckle Knee


Size	Code	Price
SMALL	MR20760AIRS	\$ 157.50
LARGE	MR20760AIRL	\$ 191.00



Restair Knee Orthosis



SIZE	TTK	KTC	THIGH CIRC	CALF CIRC
Child	12.5cm	15cm	20cm - 30.5cm	18cm - 25.5cm
X-Small	15cm	15cm	20cm - 33cm	19cm - 28cm
Small	16cm	16cm	25.5cm - 38cm	20cm - 30.5cm
Medium	19cm	19cm	33cm - 43cm	25.5cm - 38cm
Large	21.5cm	21.5cm	40.5cm - 56cm	35.5cm - 46cm



NEUROFLEX RESTORATIVE KNEE

Uses: Restores flexion lost range of motion with prolonged low load passive stretch, with flex properties to relax joint.

Features: Adjustable dials lock in to apply mild stretch at knee joint to achieve a relaxation. Side bars flex when patient draws inward, then brings the joint back to preset position. Gradually reset dials to work joint toward normal alignment and re-lengthen shortened tissue.

NEUROFLEX KNEE SEPARATOR

Uses: Soft splint to address the tightest closed knee joint.

Features: Add resister foam blocks to inside of device to gradually increase range of motion as patient improves. Turn device with bulk between knees to encourage correct alignment while side lying. Use behind knees to decrease tone when in bed. Ideal to increase comfort for Hospice patients.

SAFE STRAIGHT KNEE

Uses: Treats neurological tone of knee.

Features: Flexible, breathable, moisture wicking padding with flexible plastic sewn into chamber to provide comfortable but sturdy/therapeutic splint. Used as immobiliser or to decrease tone. Provides prolonged low load passive stretch. Keeps limbs straight without rigid, hard splints.

RETAIR KNEE

Uses: Treats knee flexion contractures with air technology.

Features: Two air bladders are inflated with the enclosed bulb hand pump. These snugly fit into sewn in chambers at the back of the leg. When inflated they take on the shape of the leg to disperse any potential pressure, and work with the hinged double uprights to increase extension/range of motion. You control the degree with flexion and extension stops on the hinges.

Sizes: ADULT Small, Medium, Large

TURNBUCKLE KNEE

Applies gentle, low-load force to knee contractures from 45° flexion to 10° hyperextension. Turnbuckle mechanism provides excellent control, allowing the therapist to adjust force in small increments

RETAIR HIP ORTHOSIS

Uses: Comfortable and safely treats hip adduction & scissoring.

Features: Two bladders run the length of each inner thigh for a total of four, securely sewn into chambers - as they are gradually inflated with a bulb hand pump, they increase the distance between the knees and help control knee scissoring. Any potential pressure is totally dispersed over the entire length of the thighs. Plastic cuffs with 3 different lengths of attachments quickly Velcro to the inside of device to further increase hip range of motion.

Sizes: Small size fits thigh circumferences from 23 to 39cm.
Large size fits thigh circumferences from 40 to 60cm.

Size	Left	Right	Price
X SMALL	MR20736XSL	MR20736XSR	\$ 56.25
SMALL	MR20736SL	MR20736SR	\$ 56.25
MEDIUM	MR20736ML	MR20736MR	\$ 56.25
LARGE	MR20736LL	MR20736LR	\$ 56.25



Restorative Hand

Size	Left	Right	Price
INFANT	MR20717IL	MR20717IR	\$ 78.75
TODDLER	MR20717TL	MR20717TR	\$ 78.75
CHILD	MR20717CL	MR20717CR	\$ 78.75
YOUTH	MR20717YL	MR20717YR	\$ 78.75



Hyper Hand

Size	Left	Right	Price
TODDLER	MR20720TL	MR20720TR	\$ 78.75
CHILD	MR20720CL	MR20720CR	\$ 78.75
YOUTH	MR20720YL	MR20720YR	\$ 78.75
ADULT ONE SIZE	MR20721L	MR20721R	\$ 78.75



Neuroflex Flex Hand

Size	Left	Right	Price
TODDLER	MR20740TL	MR20740TR	\$ 78.75
CHILD	MR20740CL	MR20740CR	\$ 78.75
YOUTH	MR20740YL	MR20740YR	\$ 78.75
MEDIUM	MR20737AML	MR20737AMR	\$ 78.75
LARGE	MR20737ALL	MR20737ALR	\$ 78.75



Thumb Ease

Size	Left	Right	Price
SMALL	MR20719SL	MR20719SR	\$ 45.00
MEDIUM	MR20719ML	MR20719MR	\$ 45.00
LARGE	MR20719LL	MR20719LR	\$ 45.00



Thumb Spica

Size	Left	Right	Price
X SMALL	MR20739XSL	MR20739XSR	\$ 56.00
SMALL	MR20739SL	MR20739SR	\$ 56.00
MEDIUM	MR20739ML	MR20739MR	\$ 56.00
LARGE	MR20739LL	MR20739LR	\$ 56.00



Cock-Up Wrist Splint

Size	Left	Right	Price
One Size	MR20715L	MR20715R	\$ 60.75



Bend Ease

Size	Left	Right	Price
SMALL	LZ52395BL	LZ52395BR	\$ 195.00
MEDIUM	LZ52395CL	LZ52395CR	\$ 195.00
LARGE	LZ52395DL	LZ52395DR	\$ 195.00



Turnbuckle Functional Position Splint

Sizing: Approximate. Small fits smaller women, Medium fits larger women and smaller men, Large fits most men



SIZE	BTE	ETF	BICEP CIRC	FOREARM CIRC
Child	9cm	10cm	16.5cm - 21.5cm	12.5cm - 18cm
Youth	10cm	12.7cm	16.5 - 21.5cm	14cm - 19cm
X-Small	10cm	12.5cm	20cm - 25.5cm	16.5cm - 21.5cm
Small	12.5cm	15cm	20cm - 30.5cm	18cm - 25.5cm
Medium	15cm	15cm	20cm - 33cm	19cm - 28cm
Large	16.5cm	16.5cm	25.5cm - 38cm	20cm - 30.5cm

NEUROFLEX RESTORATIVE ELBOW

Uses: Treats mild to severe contractures of the arm. Restore flexion lost range of motion with prolonged low load passive stretch, with flex properties to relax joint.
Features: Adjustable dials lock in to apply mild stretch at elbow joint to achieve relaxation. Side bars flex when patient draws inward, then brings the joint back to preset position. Gradually reset dials to work joint toward normal alignment and re-lengthen shortened tissue.

NEUROFLEX ELBOW SEPARATOR

Uses: Soft splint to address the tightest closed elbow joint.
Features: Add resister foam blocks to inside of device to gradually increase range of motion as patient improves. Turn device with bulk under arm to gradually and comfortably abduct shoulder and decrease tone. Use to maintain range when Neuroflex Elbow is not on. Ideal to increase comfort for Hospice patients.

SAFE STRAIGHT ELBOW

Uses: Treats neurological tone of elbow.
Features: Flexible, breathable, moisture wicking padding with flexible plastic sewn into chamber to provide comfortable but sturdy/therapeutic splint. Used as immobiliser or to decrease tone. Provides prolonged low load passive stretch. Keeps limbs straight without rigid, hard splints.

TURNBUCKLE ELBOW

Applies dynamic low-loaded force to increase elbow extension or flexion. Controls elbow contracture in patients with paralysis or significantly unbalanced muscle strength, also good for elbow flexion and extension contracture postoperatively.

HENSINGER HEAD SUPPORT


Helps children to maintain mid-line head positioning and also improve breathing, swallowing, feeding and social interaction. Must be used in conjunction with an anterior support.

HEADMASTER COLLAR

Provides excellent support and control, open design is cool and unobtrusive to wear. Bend manually for custom fit, soft rear strap with hook and loop closure.

KENTUCKY COLLAR

Uses: Effective for conditions like Torticollis/wry neck.
Features: Comfortably, safely correct side leaning of head/neck to the left. Design allows tension to relax. Remould to gradually correct neck posture.



Size	Code	Price
PAEDIATRIC	MR20724C	\$ 139.50
X SMALL	MR20724XS	\$ 139.50
SMALL	MR20724S	\$ 139.50
MEDIUM	MR20724M	\$ 139.50
LARGE	MR20724L	\$ 139.50



Neuroflex Restorative Elbow

Code	Price
MR20724SP	\$ 23.00



Neuroflex Elbow Separator

Size	Length	Circ	Code	Price
SMALL	280mm	180mm	MR20741S	\$ 45.00
MEDIUM	280mm	305mm	MR20741M	\$ 45.00
LARGE	380mm	305mm	MR20741L	\$ 45.00



Safe Straight Elbow

Overall Length	Forearm Circ	Bicep Circ	Code	Price
28.5cm	20cm - 25cm	22.8cm - 28cm	LZ52396A	\$ 195.00
36.8cm	20cm - 25cm	22.8cm - 28cm	LZ52396	\$ 195.00
36.8cm	25cm - 31.7cm	28cm - 34cm	LZ52396	\$ 195.00
41.2cm	25cm - 31.7cm	28cm - 34cm	LZ52396	\$ 195.00
41.2cm	31.7cm - 38.1cm	34cm - 40.6cm	LZ52396	\$ 195.00



Turnbuckle Elbow

Size	Circ	Code	Price
EXTRA SMALL	23cm	SPC9245-50	\$ 299.00
SMALL	30.5cm	SPC9245-51	\$ 299.00
MEDIUM	37cm	SPC9245-52	\$ 299.00
LARGE	40.5cm	SPC9245-53	\$ 299.00



Hensinger Head Support

Size	Chin to corner of jaw	Code	Price
INFANT	6.25cm	LZ52368I	\$ 295.00
CHILD	7.5cm	LZ52368C	\$ 295.00
YOUTH	8.75cm	LZ52368Y	\$ 295.00
SMALL	10cm	LZ52368S	\$ 295.00
MEDIUM	11.25cm	LZ52368M	\$ 295.00
LARGE	12.5cm	LZ52368L	\$ 295.00



Headmaster Collar

Size	Forehead to Shoulder	Code	Price
NEWBORN	> 7.5cm	MR20735NB	\$ 124.00
INFANT	7.5cm - 10cm	MR20735I	\$ 124.00
CHILD	10cm - 12.5cm	MR20735C	\$ 124.00
YOUTH	12.5cm - 15cm	MR20735Y	\$ 124.00
ADULT Large	15cm - 23cm	MR20735RAL	\$ 124.00
ADULT X Large	< 23cm	MR20735RAXL	\$ 124.00



Kentucky Collar

Please choose size and Left or Right Support.