

## KIDSTAND III / GRANSTAND III SCRIPT FORM

Customer Name: \_\_\_\_\_

Account Number: \_\_\_\_\_

Delivery Address: \_\_\_\_\_

PO Number: \_\_\_\_\_

Contact Person: \_\_\_\_\_

Contact Number: \_\_\_\_\_



### Options & Accessories

A range of accessories allows the Therapist to configure the Granstand to the individual child's clinical requirements.



Standard Foot Plates (Not Mobile)



Height Adjustable Foot System



Multi-Adjustable Foot System



Quad Pump Handle



Multi-Adjustable Knee System



Roho Knee Pads



Steel Caster Upgrade



Two Piece Handle



Contoured Body Pad with Lateral Supports and Quick Release Padded Support Belt.



Contoured Body Pad



Contoured Body Pad with Lateral Supports



Product Sizing	Kidstand	Granstand
Age — Years approx	5-15	13 onward
Max User Weight	113kg	113kg
User Height	91-142cm	157-186cm
Base Width	70cm	70cm
Base Length	102cm	102cm



# Order Form

ARTG: 163734

ITEM	ORDER CODE	QTY	PRICE
<b>Granstand III Modular Standing System</b> (Height Range: 167 - 186cm)	420C121=00000		\$ 4,500.00
<b>Kidstand III Modular Standing System</b> (Height Range: 106 - 142cm)	420C122=30000		\$ 4,300.00
<b>Pump Handle</b>			
Left			INCL
Right			INCL
<b>Frame Colour</b>			
Midnight Purple			INCL
<b>Electric Stand</b>			
Mobile Kit Granstand III: Requires Mobile Foot System	421C121=ST021		\$ 1,300.00
Electric Lift (Only on Granstand)			\$ 4,057.00
Mobile Kit Kidstand III: Requires Mobile Foot System	421C122=ST053		\$ 1,300.00
<b>Seat Sling Options</b>			
Sheepskin Sling Cover	421C121=ST032		\$ 145.00
Pocket Seat Sling Cover: Allows for additional padding of the standard sling	421C121=ST028		\$ 210.00
Vinyl Seat Sling	421C121=ST037		\$ 240.00
Contoured Vinyl Sheepskin Lined Seat Sling	421C121=ST038		\$ 290.00
<b>Lower Extremity Positioning Options</b>			
Standard Foot Pads: Width and Depth adjustable (Products are not mobile in use with this term)			No Charge
Height Adjustable Mobile Foot System: Adjust in width, height, depth, toe in and toe out. Matches the footplate height of the WC footplate. Supplies maximum growth potential, allows the Granstand III to be mobile while in use.	421C121=ST020		\$ 545.00
Multi-Adjustable Knee System: Tool-free mechanism easily adjusts in height, width and depth for optimum lower body positioning. Independent width positioning allows for proper placement of the knee pas and fast adjustment for multiple client use.	421C121=ST024		\$ 530.00
Multi-Adjustable Mobile Foot System: Individual plantar-flexion and dorsi-flexion adjustable.	421C121=ST022		\$ 725.00
Roho Knee Pads	421C121=ST030		\$ 505.00
Foot Straps	421C121=ST019		\$ 65.00
D-Ring Foot Straps	421C121=ST016		\$ 95.00
Foot Risers	421C121=ST018		\$ 187.00
<b>Upper Positioning Options</b>			
Padded Chest Belt	421C121=ST026		\$ 187.00
Small Body Pad (10.1 x 15.2cm)			\$ 56.00
Contoured Body Pad	421C121=ST013		\$ 240.00
Contoured Body Pad with Lateral Supports	421C121=ST014		\$ 500.00
Contoured Body Pad with Lateral Supports and Padded Support Belt			\$ 572.00
<b>Base/Frame Options</b>			
Steel Caster Upgrade (For hydraulic version only)	421C121=ST011		\$ 254.00
Quad Pump Handle (For hydraulic version only)	421C121=ST029		\$ 187.00
7.6cm Wider Frame			\$ 862.00
Two Piece Pump Handle			\$ 145.00
<b>Tray Options</b>			
Pneumatic Tray Lift	421C121=ST027		\$ 326.00
Tray Easel	421C122=ST056		\$ 330.00
Tray Edging	421C121=ST035		\$ 138.00
Opaque Tray Overlay	421C121=ST025		\$ 187.00
Clear Tray Overlay	421C121=ST012		\$ 187.00
Padded Tray Cover			\$ 288.00