

KIDWALK I & II SCRIPT FORM

Customer Name: _____

Account Number: _____

Delivery Address: _____

PO Number: _____

Contact Person: _____

Contact Number: _____



Standard Features

Auto Centering Pelvic Guides
 Lateral Weight Shift
 Dynamic Weight Relief
 Angle Adjustable Mast

Anti Tipper
 Wheel Locks
 Quick Release Wheels
 Padded Seat

Adjustable Axle Position
 Light Up Rear Casters
 Two Piece Frame for Easy Transport
 Tools Free Adjustability

Options & Accessories

A range of accessories allows the Therapist to configure the KidWalk to the individual child's clinical requirements.



Swing Open Lateral Supports



Upper Body Soft Support Wrap



Head Pad



1/2" Square Tube Mount



Contoured Headrests



Extended Hip Pads



Cruiser Wheels



Dynamic Leg Wraps



Extended Seats



Load Compensating Adjustment Spring



Product Sizing	KidWalk I	KidWalk II
Max User Weight	30kg	50kg
Seat Height-Inside Leg	28-55cm	37.5-70cm
Base Width	60cm	68.8cm
Base Length	67.5cm	80cm
Mast Height	82.5cm	93cm
Wheels	50cm	60cm
Castors	10cm	12.5cm



ITEM	KIDWALK I ORDER CODE	QTY	PRICE
KIDWALK I - Standard Size (W 13.9 - 22.8cm D 8.8 - 15.2cm)	420R119=00000		\$ 4,810.00
KIDWALK II - Standard Size (W 17.7 - 30.4cm D 12.7 - 17.7cm)	420R119=10000		\$ 5,400.00
Frame Colour			
Razzberry			INCL
Shocker Green			INCL
Upper Body Support			
Kidwalk I - Lateral Chest System (W 16.4 - 26.6cm D 10.1 - 16.4cm)	421R119=ST077		No charge
Kidwalk II - Lateral Chest System (W 21.5 - 36.7cm D 13.9 - 21.5cm)	421R119=ST093		No charge
<i>Includes: Lateral Side Supports, Swing away Anterior Support Pads, Lumbar Support Pad, Flip up Pad Brackets Straight</i>			
Upper Body Soft Support Wrap	421R119=ST091		\$ 325.00
Headrest Options			
Head Pad 7.6 X 10cm Flat	421R119=ST068		\$ 255.00
Contoured Headrest - Small (10cm H x 15cm W)	421R119=ST085		\$ 365.00
Contoured Headrest - Medium (12.7cm H x 20.5cm W)	421R119=ST080		\$ 449.00
1.2 cm Square Tube Headrest Mount (Headrest not included)	421R119=ST088		\$ 124.00
i2i Flip Down Headrest Mount (Headrest not included)	421C124=ST140		\$ 455.00
Lower Extremity Positioning System			
Dynamic Leg Wraps (Assist in supplying feedback to the client while encouraging movement in the lower extremities)	421R119=ST078		\$ 163.00
Extended Hip Pads - Standard on 12.7cm - 17.7cm depth	421R119=ST061		\$ 124.00
Sliding Ankle Prompt			\$ 400.00
Seat Options			
Wide Seat 7.5 cm X 20.3 cm Large Padded Seat	421R119=ST092		\$ 228.00
Leg Straps - Restricts movement of the lower extremities (Requires 15.2cm or 20.3cm extended seat)			\$ 228.00
Additional Standard Seat (ordered with unit)	421R119=ST075		\$ 130.00
Extended Seats - Longer Sides Assist in Controlling Mild Adduction and Scissoring			
Extended Seat Short (10.1cm Length)	421R119=ST065		\$ 358.00
Extended Seat Medium (15.2cm Length)	421R119=ST066		\$ 358.00
Extended Seat Long (20.3cm Length)	421R119=ST067		\$ 358.00
Frame Options			
Cruiser Wheels (Allows use as an Omni directional Walker - available for Kidwalk I Only)	421R119=ST063		\$ 389.00
Load Compensating Adjustment Spring - Supplies adjustment for bounce and client weight			\$ 384.00
Kidwalk I Wheel Options			
Spoke Wheel 50.8cm	421R119=ST059		\$ 525.00
Spoke Protectors 50.8cm	421R119=ST086		\$ 224.00
Kidwalk II Wheel Options			
Anti Roll Back Device	421R119=ST062		\$ 384.00
Spoke Wheel 60.9cm	421R119=ST060		\$ 525.00
Spoke Protectors 60.9cm	421R119=ST087		\$ 224.00
Kidwalk Directional Locks	421R119=ST095		\$ 340.00