

# SUPERSTAND

Pediatric Multi-Position Standing System

# SUPERSTAND YOUTH

Client: \_\_\_\_\_ Age: \_\_\_\_\_  
 Height: \_\_\_\_\_ Weight: \_\_\_\_\_ Dx: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Facility: \_\_\_\_\_  
 Therapist: \_\_\_\_\_ Phone: \_\_\_\_\_

PO: \_\_\_\_\_ Date: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, St, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 E-mail: \_\_\_\_\_

## SUPERSTAND

Multi-Position Pediatric Standing System

## SUPERSTAND YOUTH

Standard Features: Prone, Vertical and Supine capable; Hip and chest positioning pads with laterals; Headrest; Chinrest; Cotton support harness; Hands free standing angle adjustment; Width, height and depth adjustable knee system; Soft hip support pad; Individual planter-flexion, dorsi-flexion, rotation and height adjustable foot system; Adjustable heel loops; Toe straps; Adjustment lock out block and Base with 7.6cm plastic casters.

<input type="checkbox"/> <b>SUPERSTAND</b> Client Height: Supine <b>65cm - 120cm</b> Prone <b>75cm - 125cm</b> <b>\$3,510.00</b> Weight Limit 57 kgs Base: 52.5cm wide x 100cm deep		<input type="checkbox"/> <b>SUPERSTAND YOUTH</b> Client Height: <b>100 - 450cm</b> <b>\$4,550.00</b> Weight Limit 68 kgs Base: 60cm wide x 120cm deep
--	---	---

<b>Frame Color:</b> <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Midnight Purple <b>Harness Color:</b> <input type="checkbox"/> Black Denim <input type="checkbox"/> Blue Denim <input type="checkbox"/> Black Vinyl <input type="checkbox"/> Blue Vinyl <input type="checkbox"/> Red Vinyl	<b>Upholstery Color:</b> <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Black <b>Custom Upholstery:</b> <input type="checkbox"/> Customer Supplied Material <input type="checkbox"/> Requested Color _____
--	---

**Lateral Pads and Brackets: (All measurements in centimetres)**

Superstand is standard with 8.8H X 8.8D Chest Laterals and 13.9H X 10D Hip Laterals with 90 Degree Brackets.  
 Superstand Youth is standard with 16.4H X 11.3D Chest Laterals and 16.4H X 11.3D Hip Laterals with 90 Degree Brackets.

<b>Chest Pad Brackets:</b> <input type="checkbox"/> 90° S td. <input type="checkbox"/> 2.5cm Offset <input type="checkbox"/> 5cm Offset Upgrade <b>Hip Pad Brackets:</b> <input type="checkbox"/> 90° Std. <input type="checkbox"/> 2.5cm Offset <input type="checkbox"/> 5cm Offset Upgrade
---

<b>Superstand</b> 90 Degree:    Width adjustment from 23-33cm 2.5cm Offset: Width adjustment from 18-28cm 5cm Offset:    Width adjustment from 13-23cm	<b>Superstand Youth</b> 90 Degree:    Width adjustment from 28-38cm 2.5cm Offset: Width adjustment from 23-33cm 5cm Offset:    Width adjustment from 18-28cm
---	---

Options	<input type="checkbox"/> <b>Upper Extremity Positioning Tray Assembly:</b> Multi-Position tray with individual height, depth and angle adjustment. ....\$ <b>630.00</b>
	<input type="checkbox"/> <b>Tray Options:</b> <input type="checkbox"/> Opaque tray overlay.....\$ <b>195.00</b> <input type="checkbox"/> Clear tray overlay.....\$ <b>180.00</b> <input type="checkbox"/> Tray easel.....\$ <b>325.00</b>
	<input type="checkbox"/> <b>Activity Bowl:</b> Tray cut-out with 18cm bowl.....\$ <b>180.00</b>
	<input type="checkbox"/> <b>Padded Tray Cover:</b> Removable padded cover.....\$ <b>290.00</b>
	<input type="checkbox"/> <b>Size-Adjustable Knee Pads (pair):</b> Allows knee pad width adjustment (from 6cm - 11cm) for growth.....\$ <b>INCL</b>
	<input type="checkbox"/> <b>D-Ring Footstraps:</b> .....\$ <b>95.00</b>
	<input type="checkbox"/> <b>Shoe Holders (pair):</b> Formed ABS shoe holders with ankle and toe straps for foot positioning.....\$ <b>233.00</b> Size: <input type="checkbox"/> Sm (15cm long x 11cm wide) <input type="checkbox"/> Med (20cm long x 12cm wide) <input type="checkbox"/> Lrg (26cm long x 15cm wide)
	<input type="checkbox"/> <b>Height Adjustable Pommel:</b> Assists in maintaining lower extremity position.....\$ <b>248.00</b>
	<input type="checkbox"/> <b>Pommel W/ Rigid Hip Pad:</b> Height adjustable pommel with height and depth adjustable rigid hip pad.....\$ <b>530.00</b>
	<input type="checkbox"/> <b>Steel Caster Upgrade:</b> High performance easy rolling steel casters with lock for both swivel and roll.....\$ <b>250.00</b>
<input type="checkbox"/> <b>Padded Head Strap:</b> Sheepskin-lined strap with Velcro closure.....\$ <b>111.00</b>	



**Custom Support Pads:** Use this area only if you need custom-sized support pads on your Superstand. All custom pads must be quoted for price and manufacturing time.

Chest Pad: \_\_\_cm Wide x \_\_\_cm High Standard Pad is 25cmH x 20cmH for Superstand and 30cmW x 38cmH for Superstand Youth.


Hip Pad: \_\_\_cm Wide x \_\_\_cm High Standard Pad is 25cmW x 18cmH for Superstand and 30cmW x 20cmH for Superstand Youth.

Chest Laterals \_\_\_cm Wide x \_\_\_cm High Standard Pad is 9cmD x 9cmH for Superstand and 11cmD x 16cmH for Superstand Youth.


Hip Laterals \_\_\_cm Wide x \_\_\_cm High Standard Pad is 10cmD x 14cmH for Superstand and 11cmD x 16cmH for Superstand Youth.

**Custom Tray Surface:**

Tray: \_\_\_cm Wide x \_\_\_cm Deep Standard Tray is 57cmW x 48cmD for Superstand and 71cmW x 61cmD for Superstand Youth



Please supply any information for cutouts or modifications desired in this area or submit a paper or cardboard template. Additional information may be added in the special needs area below.



**Client Measurements:**

A: Thorax Width \_\_\_\_\_ cm      B: Hip Width \_\_\_\_\_ cm      C: Hip-to-Knee \_\_\_\_\_ cm  
 D: Knee-to-Heel \_\_\_\_\_ cm      E: Shoe Size \_\_\_\_\_      F: Shoulder Height \_\_\_\_\_ cm

Note: All Measurements Assume Symmetrical Standing and Sitting Posture.

**Special Needs:**

We are proud of our proven commitment to serving the needs of our customers with the finest standing products available since 1984. We encourage you to call us with your special needs, questions, or any item that we may assist with in your standing program. Our standing and positioning specialists will be happy to answer your questions or assist you with the design and fabrication of a unit that meets your needs. Please supply as much of the information in the client measurement area of this form as possible and list special needs in the following area to speed up the design process and / or response time.

---



---



---



---



---