

## SYMMETRY/SYMMETRY YOUTH SCRIPT FORM

Customer Name: \_\_\_\_\_

Account Number: \_\_\_\_\_

Delivery Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

PO Number: \_\_\_\_\_

Contact Person: \_\_\_\_\_

Contact Number: \_\_\_\_\_



### Standard Features

- Padded Sliding Seat
- Flip Up Depth Adjustable Kneepad
- Hydraulic Pump
- Adjustable Height Toe-In; Toe Out; Plantar-Flexion/Dorsi-Flexion Shoe Holders
- Adjustable Chest Pad
- Height and Depth Adjustable Table
- Flip Back Arm Rests
- Transport Base with 3" Plastic Casters

### Options & Accessories

A range of accessories allows the Therapist to configure the Symmetry to the individual child's clinical requirements.



Low Back



Extended Back



Lateral Supports



Swivel Seats



Hip Guides



Standard Kneepads



Headrest



Mobile Kit



Support Vest



Independent Adjustable Knee Pads



Steel Caster Upgrade



Symmetry shown with optional low back, extended back, hip guides, lateral supports, two piece handle and headrest.



Solid Calf Panel



Roho Knee Pads



Product Sizing	Symmetry Youth	Symmetry
Age — Years approx	5-15	13 onward
Max User Weight	113kg	124.5kg
User Height	106-167cm	152-195cm



# Order Form

ARTG: 163732

ITEM	ORDER CODE	QTY	PRICE
<b>Symmetry</b> (User Height: 152.4 - 194.7cm User Weight Limit: 124.5kg)	420C124=80000		\$ 4,500.00
<b>Symmetry Youth</b> (User Height: 106.6 - 167.6cm User Weight Limit: 113.2kg)	420C124=90000		\$ 4,500.00
<b>Pump Handle</b>			
Left			INCL
Right			INCL
<b>Upholstery Colour</b>			
Grey			INCL
Blue			INCL
<b>Seat Options</b>			
Swivel Seat (Seat swivels and locks in multiple positions)	421C124=ST136		\$ 860.70
Hip Guides (Provides hip and lower extremity positioning and alignment)	421C124=ST124		\$ 347.00
Seat Belt (Velcro D-Ring Belt)	421C124=ST133		\$ 72.50
<b>Back Options</b>			
Back systems feature our adjustable sliding back system. Our backs move with the client eliminating the effect of shear. Each back adjusts allowing fine tuning of the back support items for the individual user when desired.			
Low Back 22.8cm tall	421C124=ST127		\$ 475.00
Extended Back; Extension mounts to low back to create a 50.8cm high back	421C124=ST120		\$ 260.00
Lateral Supports: 15.2cm x 15.2cm triangle shaped pads for width and height and move with the sliding back (requires back assembly)	421C124=ST126		\$ 360.00
Headrest: Height and Depth adjustable for support while sitting or standing	421C120=ST140		\$ 355.00
Support Vest: Anterior support vest in denim	421C123=ST109		\$ 280.00
Chest Belt: Single strap belt attaches to frame with Velcro loop	421C124=ST118		\$ 115.00
Padded Chest Belt: Single strap with movable pad attached to frame with Velcro loop	421C124=ST130		\$ 188.00
<b>Lower Extremity Positioning Options</b>			
Independent Adjustable Knee Supports: Independent depth adjustable allows for additional lower extremity positioning	421C124=ST125		\$ 460.00
Roho Knee Pads			\$ 500.00
Foot Straps: Velcro straps to assist foot positioning	421C124=ST122		\$ 65.00
D-Ring Foot Straps: Velcro straps with D-Rings for additional holding strength			\$ 94.00
Solid Calf Panel: Height adjustable padded panel attaches under the seat to maintain the feet and lower extremities in a neutral position	421C124=ST117		\$ 200.00
<b>Frame Options</b>			
Steel Caster Upgrade	421C123=ST108		\$ 252.00
Mobile Kit	421C124=ST128		\$ 1,300.00
Two Piece Pump Handle	421C123=ST115		\$ 145.00
Quad Pump Handle	421C124=ST132		\$ 145.00
Tray Edging	421C123=ST114		\$ 136.00
Tray Easel: Adjustable angle easel sits on tray surface	421C123=ST113		\$ 315.00
Padded Tray Cover	421C124=ST131		\$ 290.00
Clear Tray Overlay	421C124=ST119		\$ 188.00
Opaque Tray Overlay	421C124=ST129		\$ 188.00