

Granstand

Modular Standing System

One Stander for Life!

Standing you the way
nature intended.



Features

- No transfer, no lifting design
- Independent mobility
- Growth capable
- Space saving expandable frame
- Advanced foot positioning and mobility
- Upgradeable modular frame construction
- Hydraulic or Electric Standing system - You choose
- Kidstand only available in Hydraulic



Kidstand

Granstand

Complete Hip Extension, Contracture Control, Progressive Standing



Head Office: 2/19 Boden Road, Seven Hills, NSW 2147
Phone: 02 9674 8904
Fax: 02 8212 5912
Website: www.dmedirect.com.au

Granstand

Modular Standing System

ARTG 163734

Product Sizing	Kidstand	Granstand
Age — Years approx	5-15	13 onward
Max User Weight	113kg	113kg
User Height	91-142cm	157-186cm
Base Width	70cm	70cm
Base Length	102cm	102cm



Mobile System

The Granstand aligns the body in its natural standing position with a totally individualised fit. This is accomplished through the most extensive system of adjustable components available. Individual chest pad, tabletop, kneepad and footpad adjustments along with various sling options assure the best standing position possible.

The innovative design features of the Granstand allow the necessary adjustments to facilitate complete hip extension. This comfortable, personalised fit reduces stress at pressure points and supplies proper weight bearing posture providing maximum weight bearing as well as contracture control.

Our modular system allows for unsurpassed growth ability at a minimum cost both in time and money. Most growth adjustments require only a few minutes and many of these changes require the purchase of little or no additional equipment. If the clients growth is accompanied by a change in physical condition, equipment upgrades can be added as needed without purchasing a new standard frame.



Expanding Leg System With Standard Leg Rest

Options & Accessories

A range of accessories allows the Therapist to configure the Granstand to the individual client's clinical requirements.



Standard Foot Plates (Not Mobile)



Height Adjustable Foot System



Multi-Adjustable Foot System



Multi-Adjustable Knee System



Contoured Body Pad with Lateral Supports and Quick Release Padded Support Belt.



Steel Caster Upgrade



Contoured Body Pad



Contoured Body Pad with Lateral Supports



Two Piece Handle



Roho Knee Pads



Quad Pump Handle