

KidWalk

Dynamic Mobility System

Because Life is Full
Of Important Steps!



Features

- Quick exacting pelvic positioning
- Single adjustment bilateral pelvic pad placement
- Mid wheel location for better control and smaller turning radius
- Swing up lateral supports
- Swing up pelvic supports
- Hip and chest laterals adjustable in height, width and depth



A New Dimension in Mobility...
Designed for clients to access and explore their environment!



Durable Medical Equipment Pty. Ltd.



Head Office: 2/19 Boden Road, Seven Hills, NSW 2147
Phone: 02 9674 8904
Fax: 02 8212 5912
Website: www.dmedirect.com.au

KidWalk


Dynamic Mobility System

ARTG 162286

Product Sizing	KidWalk I	KidWalk II
Age — Years approx	1-7	5-12
Max User Weight	30kg	50kg
Seat Height-Inside Leg (1)	28-55cm	37.5-70cm
Hip Width (2)	14-23.5cm	18-29cm
Chest Width (3)	16.5-25.5cm	22-37cm
Thoracic Height (4)	24.5-41cm	30-56cm
Base Width	60cm	68.8cm
Base Length	67.5cm	80cm
Mast Height	82.5cm	93cm
Overall Length	82.5-87.5cm	90-95cm
Wheels	50cm	60cm
Castors	10cm	12.5cm



Dynamic vertical movement allows for stepping phase



Dynamic lateral movement allows for user weight shift



KidWalk disassembled for transport


The KidWalk is the first partial weight bearing, posturally supportive mobility system that is designed to move with the child, giving them the freedom to walk unaided while keeping their hands free to touch, feel and learn from their surroundings.

KidWalk Suitable for mild, moderate and complex users with lower limb motor function. A weight shift mechanism allows work on balance and motor development during therapy and a mid-wheel design encourages rotation of the upper body over the pelvis while turning.

The small turning circle and ability to roll over carpet make it suitable for inside the home as well as the classroom, therapy room and outdoors.

Options & Accessories

A range of accessories allows the Therapist to configure the KidWalk to the individual child's clinical requirements.



Wheel Locks



Cruiser Wheels



Hip Laterals



Saddle Options



Head Pad



Head Rest Options



Laterals



Seat Mounted Leg Straps



Leg Wraps



Wheels



Head Rest Mounts



Chest Strap