



# DME

*Durable Medical Equipment Pty. Ltd.*


## Product Catalogue

# Spare Parts





## ARM REST PADS



### Feature

- 6 Fixing points
- High density foam
- Plastic Shell
- Lightweight - Desk 150g, Full 232g

### Benefit

- Will fit most armrests
- Hard wearing
- Clean Lines/ Professional finish
- Doesn't add unnecessary weight to mobility system

Code	Description	
MIAAR005-HR	Desk Arm Rest Pad	EACH
MIAAR003-HR	Full Arm Rest Pad	EACH

## SEAT BELTS (Positioning)



### Feature

- 50mm Black Webbing
- Plastic or Auto Buckle
- Fixing Toggles (wheelchair belt)
- 14mm Fixing Grommet (power chair/scooter belt)
- Wheelchair belt adjusts to 1.52m
- Power chair/Scooter belt adjusts to 1.2m

### Benefit

- Strong and Durable
- Suit users needs
- Easy to fit to any style of chair
- Secure strong attachment
- Fits most users

Code	Description
MIASB-DMEP01	Wheelchair - plastic buckle
MIASB-DMEA01	Wheelchair - Auto Buckle
MIASB001	Power Chair/ Scooter Plastic Buckle
MIASB109	Power Chair/ Scooter Auto Buckle



## UPHOLSTERY

### Feature

Manufactured using Black 600 Polyester  
 Back Upholstery features back pocket 250x200mm with flap  
 Rip Stop Inner

**NB: Upholstery NOT supplied with rails or fixing screws**

### Benefit

Strong, durable  
 Useful for storage  
 Reduces Wear & Tear  
 Easy to Clean



Code	Description
MIABS-ST1616	SEAT 41 x 41cm
MIABS-ST1618	SEAT 41 x 46cm
MIABS-ST1816	SEAT 46 x 41cm
MIABS-ST1818	SEAT 46 x 46cm
MIAS-BK1617	BACK 41 x 43cm
MIAS-BK1817	BACK 46 x 43cm

## BRAKES


### Feature

U style Clamp fits 22mm tube  
 Push to Lock  
 Large Brake Shoe  
 Chrome Steel construction

**Sold in pairs**

### Benefit

Fits Most Standard Chairs  
 Easy for Users to Operate  
 Sits firmly on Wheel  
 Strong and Durable



Code	Description	
MHP8665/1611A	Brake Self Propelling 7cm brake arm	PAIR
MHP8665/1611B	Brake Transit 10cm brake arm	PAIR

## HAND GRIPS

### Feature

Ergonomic Shape  
 Fits 22m tube

### Benefit

Fits Hand Snugly  
 Fits Most Standard Chairs



Code	Description	
MHP8665/2341	Hand Grip	EACH



## DME Tray Clamps

Therafin produce some of the most effective multi-fit tray clamps around!

### STANDARD ARM or TUBE ARM SLIDE ON E-Z LOCK CLAMPS



- Can be mounted flush to the outside of the armpads when tray width is an issue
- Durable glass-filled nylon black material features recessed hardware that won't tear armpads
- Slide the tray on from the front
- Clamp-on design securely holds tray
- Armpad height and width adjustable
- Sold in pairs, attaching hardware for 0.95cm tray included



### TOP DROP CLAMPS



- Requiring only 0.95cm under-armpad clearance, it's ideal for use when a control system is mounted close to the armpad
- Position the tray straight down from the top without disturbing the joystick or control system
- Clamp-on design securely holds tray
- Armpad height and width adjustable
- Durable glass-filled nylon black material features recessed hardware that won't tear armpads
- Sold in pairs, attaching hardware for 0.95cm tray included

