

Roho® AirLITE® Cushion



Comfort

The Roho® AirLITE® Cushion is unlike any Cushion you've ever seen.

Designed for comfort and peace of mind, the AirLITE combines the benefits of a sealed air support pad and a durable contoured foam to provide postural support and optimal durability.

Positioning

The non-adjustable pre-set air support pad, uses patented Roho Air Floatation technology, while a contoured high density foam base positions the body in a comfortable position, creating one of the most effective and easy to use cushions on the market today.

At only 0.5kg it continues the tradition of lightweight products from the Roho group.

Simplicity

AirLITE, comfort, durability and positioning in one simple package, the solution the others can't match.

Roho® AirLITE® Cushion

By combining a contoured anatomical shaped foam with a pre-set, sealed air floatation component the Roho® AirLITE® Cushion enables clients to sit comfortably for an extended period of time. AirLITE's technology and design provides the comfort clients want and the support they need.

Most foam cushions will compress over time, but the AirLITE's pre-set, sealed Air Floatation component provides comfort and peace of mind you can count on long after other "just foam" cushions have lost their effectiveness.

The contoured shape of the AirLITE cushion provides clients with support for basic positioning needs. Leg troughs and an ischial well provide support for lower extremity control, stability and comfort.

The AirLITE cushion weights only .98kg. Compare this to alternative technologies at two or four times the weight of the AirLITE. A lightweight construction makes a mobility system easier to propel and transport.



FEATURES AND BENEFITS

Simplicity: It is easy - No adjustment necessary

Washable: Dry Floatation support pad can be easily cleaned with soap and water

Security: The non-skid bottom increases security of the cushion in the chair

Foam Base

Fluid Resistant: Breathable, fluid resistant cover protects the cushion from odour and is machine washable

SPECIFICATIONS

Includes : Airlite Cushion, 2 way stretch cover, Operating Instructions

Construction: Flame resistant polyurethane seating foam. Pre-set air cell: polyurethane

Sizes: Range from 35cm up to 50cm

User Weight Limit: 114kg

MORE PRODUCTS



ROHO® Sofflex® 2
Mattress Overlay System



ROHO® Mosaic® Cushion