

Roho[®] High Profile[®] Cushion Dual Compartment Cushion

10cm cell height offers maximum therapeutic performance for the treatment and prevention of ischemic ulcers


Individual air cells allow air circulation between the cells

Lightweight only 1.5kgs

NO User Weight Limit

Shapes to the individual addressing their specific seating needs.

Durable Neoprene, Flame resistant and resistant to urine, oils and medicines



Durable

The Roho[®] High Profile[®] Dual Compartment Cushion offers maximum therapeutic performance for the treatment and prevention of ischemic ulcers.

Adjustable

The High Profile Dual Compartment Cushion adjusts to the shape of the individual.

Lightweight

The cushion is only 1.5kg making it easy to transport and handle, and minimises fatigue of the individual.

Therapeutic



Head Office: 2/19 Boden Road, Seven Hills, NSW 2147
Phone: 02 9674 8904
Fax: 02 8212 5912
Website: www.dmedirect.com.au

Roho® High Profile® Cushion

Dual Compartment Cushion

The Roho® High Profile® Dual Compartment Cushion provides an adjustable, low deformation, counter pressure environment that not only help prevent, but assists in the healing of ischemic ulcers.

The High Profile Dual Compartment Cushion has two inflation valves which independently adjust two separate compartments of the cushion to the needs of the user. The two separate compartments allow for therapeutic positioning of the individual while achieving low sitting pressures.

The fully adjustable cushion is lightweight and consists of 10cm high interconnected air cells which allows for maximum therapeutic performance and enhanced sitting times.

The High Profile Dual Compartment Cushion can position the user from side to side to adjust for a pelvic obliquity or from front to back to address sacral sitting.

FEATURES AND BENEFITS

10cm Cell Height:	Offers maximum therapeutic performance for the treatment and prevention of ischemic ulcers
Adjustable:	Shapes to the individual addressing their specific seating needs
Repairable:	Extends overall life of mattress
Lightweight:	Cushion is 1.5kg, making it easy to transport and handle
Washable, Two Way Stretch Cover:	Allows for proper immersion into the cushion and "breathes" with every shift in position
Durable Neoprene Rubber:	Flame resistant, easy to clean or steam autoclave


SPECIFICATIONS

- Material:** Neoprene
- Seat Widths:** Range from 30cm to 70cm Widths
- Weight:** 1.5kg
- Maximum User Weight:** NONE

MORE PRODUCTS



Jet Stream Pro® Back Support System
Ultra lightweight back rest



ROHO® Dry Floatation® Mattress Overlay System