


Roho[®] Low Profile[®] Single Compartment Cushion

5cm cell height offers maximum therapeutic performance for the treatment and prevention of ischemic ulcers

Individual air cells allow air circulation between the cells

Lightweight only 1.5kgs

NO User Weight Limit



Shapes to the individual addressing their specific seating needs.

Durable Neoprene, Flame resistant and resistant to urine, oils and medicines

Durable

Adjustable

Lightweight

Therapeutic

The Roho[®] Low Profile[®] Single Compartment Cushion offers maximum therapeutic performance for the treatment and prevention of ischemic ulcers.

The Low Profile Single Compartment Cushion adjusts to the shape of the individual.

Cushion is only 1.5kg making it easy to transport and handle, and minimises fatigue of the individual.



Head Office: 2/19 Boden Road, Seven Hills, NSW 2147

Phone: 02 9674 8904

Fax: 02 8212 5912

Website: www.dmedirect.com.au

Roho® Low Profile® Single Compartment Cushion

The Roho® Low Profile® Single Compartment Cushion provides an adjustable, low deformation, counter pressure environment that not only prevents, but assists in healing of ischemic ulcers.

The Low Profile Single Compartment Cushion remains the cushion of choice for individuals who suffer from any stage ischemic ulcer, are at high risk, or have a history of tissue breakdown.

The fully adjustable cushion is lightweight and consists of flexible 5cm high interconnected air cells which allow for maximum therapeutic performance.

The Low Profile Single Compartment Cushion, the ultimate in therapeutic cushioning.



FEATURES AND BENEFITS

5cm Cell Height:	Offers maximum therapeutic performance for the treatment and prevention of ischemic ulcers
Adjustable:	Shapes to the individual addressing their specific seating needs
Repairable:	Extends overall life of mattress
Lightweight:	Cushion is 1.5kg, making it easy to transport and handle
Washable, Two Way Stretch Cover:	Allows for proper immersion into the cushion and “breathes” with every shift in position
Durable Neoprene Rubber:	Flame resistant, easy to clean or steam autoclave


SPECIFICATIONS

Material:	Neoprene
Seat Widths:	Range from 30cm to 70cm Widths
Weight:	1.5kg
Maximum User Weight:	NONE

MORE PRODUCTS



Jet Stream Pro® Back Support System
Ultra lightweight back rest



ROHO® Dry Floatation® Mattress Overlay System