

# Roho<sup>®</sup> Mosaic<sup>®</sup> Cushion



## Skin Protection

The Roho<sup>®</sup> Mosaic<sup>®</sup> Cushion is especially suited for users who require a cushion for medium to high level skin protection, support and comfort.

## Simplicity

Using air flotation technology, the Mosaic cushion is not only comfortable, but the use of interconnected air cells to distribute forces evenly across the cushion allows it to be used effectively for at risk skin protection.

## Versatile

Cushion features drainage holes in base of cushion allowing cushion pad to be used in wet areas such as on shower seats and or shower stools.

# Roho<sup>®</sup> Mosaic<sup>®</sup> Cushion

The Roho<sup>®</sup> Mosaic<sup>®</sup> Cushion is especially suited for users who require a cushion for medium to high level skin protection, support and comfort.

Using air flotation technology, the Mosaic cushion is not only comfortable, but the use of interconnected air cells to distribute forces evenly across the cushion allows it to be used effectively for at risk skin protection.

Cushion features drainage holes in base of cushion allowing cushion pad to be used in wet areas such as on shower seats and or shower stools.

## FEATURES AND BENEFITS

**Lightweight:** Adds very little weight to mobility system - 0.5kg

**Washable:** Can be easily cleaned with soap and water. Cover is machine washable

**Repairable:** Polyvinyl Chloride cushion can be patched to extend use

## SPECIFICATIONS

**Includes:** Mosaic Cushion, 2 Way Stretch Cover, Operating Instructions, Pump, Repair Kit

**Construction:** Polyvinyl Chloride

**Sizes:** Range from 40cm up to 45cm

**User Weight Limit:** 114kg

## MORE PRODUCTS



**ROHO<sup>®</sup> Harmony<sup>®</sup> Cushion**



**ROHO<sup>®</sup> LTV<sup>®</sup> Cushion**