

# Roho® Quadtro Select® Mid Profile™ Cushion

7.5cm cell height offers therapeutic performance for the treatment and prevention of ischemic ulcers

Individual air cells allow air circulation between the cells

Shapes to the individual addressing their specific seating needs

Isoflo memory control allows you to fit client while seated

No user weight limit

Durable Neoprene, Flame resistant and resistant to urine, oils and medicines



## Positioning

The Roho® Quadtro Select® Mid Profile® Cushion sets the standard in the industry for overall performance in wheelchair seating.

## Skin Protection

Each and every user will benefit from a customised fit through the simple push of a knob. The revolutionary ISOFLO Memory Control offers shape fitting capabilities while the user is seated, allowing quick and easy, on-demand adjustment to maximise function.

## Stability

With the built in stability and simplicity of the Quadtro Select Mid Profile Cushion, no longer will you have to sacrifice maximum skin protection to get stability, positioning or convenience. The Quadtro Select Mid Profile Cushion does it all. And with no compromise.

## Convenience



Durable Medical Equipment Pty. Ltd.

Head Office: 2/19 Boden Road, Seven Hills, NSW 2147

Phone: 02 9674 8904


Fax: 02 8212 5912

Website: [www.dmedirect.com.au](http://www.dmedirect.com.au)

# Roho<sup>®</sup> Quadtro Select<sup>®</sup> Mid Profile<sup>™</sup> Cushion

The Roho<sup>®</sup> Quadtro Select<sup>®</sup> Mid Profile<sup>®</sup> Cushion allows for quick, simple position selection for every user while seated. By pushing the knob to the open position, you can adjust the users posture to the desired sitting position. By pushing the knob to the closed position you can maintain the desired sitting posture for maximum function.

Based on thirty years research and client experience patented Roho<sup>®</sup> Dry Floatation<sup>®</sup> technology provides unparalleled performance in facilitating blood flow by fitting matching and tracking the shape of a person.



The Quadtro Select Mid Profile Cushion offers a new level of stability never before found in an air floatation cushion. Stability is enhanced by locking air into place in each of the four quadrants through the simple use of the isoflo memory control, minimising side to side or front to back motion.

All postural adjustments are made while the user is seated, eliminating the need for transfer on and off the chair during the evaluation process. In addition the Quadtro Select Mid Profile Cushion allows the user to selectively change their position throughout the day as their particular activity demands.

## FEATURES AND BENEFITS

- Adjustable:** Fits specific seating needs and unique body shapes
- Lightweight:** Easy to handle and carry at just 1.7kgs
- Washable:** Can be easily cleaned with soap and water. Cover is machine washable
- Repairable:** Neoprene rubber cushion can be patched or repaired to extend use
- Ventilated:** Separate air cells allow air circulation to reach surface of the skin and provide a cooling environment

## SPECIFICATIONS

**Includes:** Quadtro Select Cushion, 2 way stretch cover, hand inflation pump, repair kit, operating instructions

**Construction:** Flame resistant, neoprene rubber


**Sizes:** Range from 30cm up to 55cm

**User Weight Limit:** NONE

## MORE PRODUCTS



**Jet Stream Pro<sup>®</sup> Back Support System**  
Ultra lightweight back rest



**ROHO<sup>®</sup> Dry Floatation<sup>®</sup> Mattress Overlay System**