


Roho[®] nexus Spirit[®] Cushion

Exceptional positioning of the pelvis and lower extremities for enhanced sitting posture

Individual air cells allow air circulation between the cells



Contoured foam base is designed for increased stability for transfers

Roho Cushion positioned to address critical pressure points provides you the highest level of skin protection

No user weight limit

Skin Protection

The Roho[®] nexus Spirit[®] Cushion is unlike any Roho Cushion you've ever seen before. It combines the stability of a contoured foam base with the therapeutic care for which Roho[®] Dry Floatation[®] cushions are known.

Simplicity

The contoured foam base is designed for increased stability for transfers and provides exceptional positioning of the pelvis and lower extremities for enhanced sitting posture.

Stability

The Roho Cushion is positioned to address critical pressure points provides you the highest level of skin protection from the originators of Roho Dry Floatation technology.

Life is full of challenges, sitting shouldn't be one of them.

Roho[®] nexus Spirit[®] Cushion

The Roho[®] nexus Spirit[®] Cushion provides protection to the ischia, sacrum and coccyx area through the placement of a Roho[®] Dry Floatation[®] support pad. Unparalleled in performance the Roho Cushion provides an adjustable, low friction, low shear environment that can protect the user against wounds.

The contoured foam base provides increased stability for transfers in and out of the cushion.

The contoured foam is designed to provide improved positioning of the pelvis and lower extremities, enhancing sitting posture.

FEATURES AND BENEFITS

Lightweight: Adds very little weight to mobility system - 1.4kg

Washable: Dry Floatation support pad can be easily cleaned with soap and water. Foam base may be wiped down. Fluid resistant cover is machine washable.

Repairable: Neoprene rubber support pad can be patched to extend use.

SPECIFICATIONS

Includes: Nexus Spirit Cushion, 2 way stretch cover, hand inflation pump, repair kit, Operating instructions

Construction: Flame resistant neoprene rubber support pad, Flames resistant polyurethane seating foam


Sizes: Range from 36cm up to 51cm

User Weight Limit: NONE

MORE PRODUCTS



Jet Stream Pro[®] Back Support System
Ultra lightweight back rest



ROHO[®] Dry Floatation[®] Mattress Overlay System