

Superstand

Multi-Positioning Standing System

One Frame delivers

- * Prone
- * Supine
- * Upright

No tool adjustment!



Features

- Converts from Prone, Upright to Supine with no tools
- Fully adjustable support pads, straps and harness
- Knee System adjusts for width, depth & height
- Foot system, height, width, depth angle and rotation adjustments
- Mobile Base
- Tray, height depth and angle adjustable
- Built in growth

Supporting you *when you need it, how you need it.*



Head Office: 2/19 Boden Road, Seven Hills, NSW 2147
Phone: 02 9674 8904
Fax: 02 8212 5912
Website: www.dmedirect.com.au

Superstand

Multi-Positioning Standing System

Product Sizing	Superstand	Superstand Youth
Age — Years approx	2-8	6-14
Max User Weight	57kg	68kg
User Height	65-125cm	100-150cm
Hip Width	13-33cm	18-38cm
Chest Width	13-33cm	18-38cm
Base Width	52.5cm	60cm
Base Length	100cm	120cm

Tool free adjustments allow easy changes to our body support pads, headrest, chinrest, footrest, knee supports, support straps, chest harness and standing angle. Wide range of standing angle adjustment allows for ease of client positioning, adjustment and weight loading. Soft padded positioning supports and cushioned straps, plus complete adjustability, properly position the client to maximise standing time and tolerance.

Superstand is engineered to grow with your standing needs. Superstand fits clients from 66-122cm and Superstand youth fits clients from 102-152cm allowing for maximum growth throughout each individual standing program, while maintaining full positioning capabilities. Multiple choices in pad mounting brackets allow for size and growth changes and orthopaedic positioning needs.



SUPINE



PRONE

Options & Accessories

A range of accessories allows the Therapist to configure the Superstand to the individual child's clinical requirements.



Shoe Holders




Steel Casters



Size Adjustable Knee Pads



Activity Bowl



Height Adjustable Pommel with Rigid Hip Pad

Tray
 Standard Tray is: Superstand: 51cm Wide x 48cm Deep & Superstand Youth: 71cm Wide x 61cm Deep

Custom Tray:
 Need to supply measurements for Width and Depth. Also supply any information for cutouts and modifications desired by providing a paper or cardboard template.