

# Symmetry

Supporting you  
when you need it,  
how you need it.

## Features

- Seat and Back Swivel as one unit
- Support Surfaces that Stay With You
- Full Upright Standing Posture
- Shear Reducing Sliding-Seat
- Adjustable-Sliding-Back
- Easier Transfers



**It's all about the Swivel !!!**



Head Office: 2/19 Boden Road, Seven Hills, NSW 2147  
Phone: 02 9674 8904  
Fax: 02 8212 5912  
Website: [www.dmedirect.com.au](http://www.dmedirect.com.au)

# Symmetry

ARTG 163732

## A new Spin in Standing Technology

Product Sizing	Symmetry Youth	Symmetry
Age — Years approx	5-15	13 onward
Max User Weight	113kg	124.5kg
User Height	106-167cm	152-195cm



Standard System

Flip-back armrests, flip-up lift arms, open frame design and flip-up knee support allow maximum space for ease of entry and exit with any type of transfer. Lateral or standing pivot transfers are made easier without any components in the way.

Symmetry's open frame design allows better access to the frame during transfers, even when a lifting device is needed. This open design supplies extra space for getting your chair in close or allowing attendants to work closely.

The Symmetry Sliding-Seat technology reduce the effect of shear as the body moves from sitting to standing. Our innovative Sliding-Seat assists in maintaining the position of the body while minimizing unwanted pressures and shear on the seated areas, allowing you to stay in position. Our Adjustable-Sliding-Back virtually eliminates the effects of shear and body movement when going from sitting to standing. Each back assembly is adjustable so it will move with the client, matching the needs of each individual.

Symmetry ability for the back and seat to maintain proper alignment while standing is unequalled. The design of the 'support system maintains the back in a full-contact supportive position, both while seated and, most importantly, while standing.



Low Back System



High Back System

### Options & Accessories

A range of accessories allows the Therapist to configure the Symmetry to the individual client's clinical requirements.



Low Back



Extended Back



Lateral Supports



Swivel Seats



Hip Guides



Standard Kneepads



Headrest




Symmetry shown with optional low back, extended back, hip guides, lateral supports, two piece handle and headrest.



Mobile Kit



Support Vest



Independent Adjustable Knee Pads



Steel Caster Upgrade



Solid Calf Panel



Roho Knee Pads