

LTV[®] Seat Cushion

Individual inter connected air cells to distribute body weight pressure

Padded Cover for additional comfort



Low profile cells for use on multiple seating surfaces.

Lightweight and easy to fold and transport

NO user weight limit

Comfort

You'll Love that Valve!

The secret to the extraordinary comfort of this cushion is the valve.

It lets you customise the seat to your specific comfort level. Roho shape fitting technology makes the difference.

Lightweight

The LTV[®] Seat[®] Cushion is made of individual, interconnected air cells that allow air to slowly transfer from chamber to chamber and evenly distribute body weight pressure.

Simplicity

Weighting only 0.45kg, the LTV Seat is lightweight and easy to fold. It can easily be transported from your scooter, to the car or wherever you may need a cushion for use outside of a wheelchair.

LTV[®] Seat Cushion

The LTV[®] Seat[®] Cushion from ROHO is a great little cushion that you can take with you anywhere!

This folding, lightweight cushion is extremely portable, & is so versatile it can be used in a variety of settings. It can be used in the office, at sporting events, at home, in the car, & on mobility scooters.

This cushion has dynamic air chambers that allow air to shift between cells, to better contour to your unique shape. The LTV Seat automatically adjusts to your movement, minimizing pressure & providing soothing comfort when you need it the most. The key to this cushion is the twist valve that allows you to make minute adjustments, & customize the air pressure to your exact liking.

This cushion is a great choice for anyone who experiences pain after long periods of sitting. The LTV seat has so many benefits, & it's available at a very low cost. How long are you willing to deal with pain & discomfort?

Choose the LTV Seat from ROHO with Shape Fitting Technology, & take a step toward living pain-free.

LTV Seat Cushion
LTV-xxxx

FEATURES AND BENEFITS

Lightweight: Only 0.45kg take it with you anywhere

Washable: Air support pad can be easily cleaned with soap and water

Padded Cover: Provides additional comfort

SPECIFICATIONS

Includes: LTV Cushion, Padded cover, Operating instructions


Size: 42.5cm x 47.5cm

User Weight Limit: NONE

MORE PRODUCTS



Jet Stream Pro[®] Back Support System
Ultra lightweight back rest



ROHO[®] Dry Floatation[®] Mattress Overlay System