

Hip Braces

RESTAIR HIP ORTHOSIS

MR207060air



Uses: Comfortable and safely treats hip adduction & scissoring.

Features: Two bladders run the length of each inner thigh for a total of four, securely sewn into chambers - as they are gradually inflated with a bulb hand pump, they increase the distance between the knees and help control knee scissoring. Any potential pressure is totally dispersed over the entire length of the thighs. Plastic cuffs w/ 3 different lengths of attachments quickly Velcro to the inside of device to further increase hip range of motion. Most comfortable. Safe, effective. Great for positioning.

Sizes: Small size fits thigh circumferences from 23 to 39cm.
Large size fits thigh circumferences from 40 to 60cm.