


# Neck Braces



SIZE	FOREHEAD TO SHOULDER
Newborn	6.5cm - 9cm
Infant	7.5cm - 10cm
Child	10cm - 12.5cm
Youth	12.5cm - 15cm
Adult	15cm - 23cm

Please choose size and Left or Right Support.

## KENTUCKY KOLLAR

MR20735



- Uses:** Effective for conditions like Torticollis/wry neck.
- Features:** Comfortably, safely correct side leaning of head/neck to the left. Design allows tension to relax. Remould to gradually correct neck posture.
- Sizes:** PAEDIATRIC, Infant, Child, Youth; ADULT, Universal (One Size)